

# SunQ Sil HPLC Columns

**Make your Chromatography Easy**



## SALIENT FEATURES:

- HIGH PERFORMANCE - TO - PRICE RATIO.
- MADE FROM ULTRA-PURE SILICA GEL WHICH ENSURES HIGH MECHANICAL STRENGTH WHILE MINIMIZING FINES.
- AVAILABLE IN VARIETY OF PARTICLE SIZES AND PORE SIZES FOR EFFICIENT MICRO AND SEMI-MICRO APPLICATIONS.
- AVAILABLE IN VARIETY OF CHEMICAL SURFACE BONDINGS THAT ENSURE REVERSE PHASE AND NORMAL PHASE CHROMATOGRAPHY.
- DURABLE AND TOLERANT EVEN WITH 100% AQUEOUS MOBILE PHASE.
- LAB TESTED WITH STRICT QA/QC FOR OPTIMUM PERFORMANCE.
- OFFERED BY SUNTEK SERVICES , INDIA.

For enquiries,

Suntek Services; E: [sales@suntek-services.com](mailto:sales@suntek-services.com); PHONE NUMBER: +91-22-2610 4285/ 9608, 6697 6816/ 17, 98670 26530

# SUN INTERNATIONAL HPLC COLUMNS

**High quality analytical columns manufactured in Japan.**

High quality analytical columns manufactured in Japan.

Below is the list of standard columns of 4.6mm ID of 250mm and 150mm length.

C18HS-5, 4.6mm ID X 250mm L, 100Å	C18HS-5, 4.6mm ID X 150mm L, 100Å
C8-5, 4.6mm ID X 250mm L, 120Å	C8-5, 4.6mm ID X 150mm L, 120Å
NH2-5, 4.6mm ID X 250mm L, 120Å	NH2-5, 4.6mm ID X 150mm L, 120Å
PH-5, 4.6mm ID X 250mm L, 120Å	PH-5, 4.6mm ID X 150mm L, 120Å
CN-5, 4.6mm ID X 250mm L, 120Å	CN-5, 4.6mm ID X 150mm L, 120Å
SIL-5, 4.6mm ID X 250mm L, 120Å	SIL-5, 4.6mm ID X 150mm L, 120Å
C4-5, 4.6mm ID X 250mm L, 120Å	C4-5, 4.6mm ID X 150mm L, 120Å
C1-5, 4.6mm ID X 250mm L, 120Å	C1-5, 4.6mm ID X 150mm L, 120Å

We at Suntek Services have every columns available with different specifications i.e 4.0mm ID or length of 35mm, 50mm, 75mm, 100mm & 300mm.

Additional we have the Cartridge Columns for the above mentioned main columns.

- > Cartridge Guard Column Adaptor.
- > Cartridge Guard Column, 4.6mm ID x 10mm L
- > Cartridge Guard Column, 4mm ID x 10mm L

We also have Preparative Columns available of every material with different specifications i.e

7.8mm, 10mm, 21.2mm, 30mm, 50mm Length: 250mm LA: 100 A

For enquiries,

Suntek Services; E: [sales@suntek-services.com](mailto:sales@suntek-services.com); PHONE NUMBER: +91-22-2610 4285/ 9608, 6697 6816/ 17, 98670 26530